

## DenApex Apex Locator



DenApex Apex Locator adopts multi-frequency operating system, which adapts into different root canal conditions. High contrast colour screen makes the display nice and clear.

### Specifications



Voltage	DC4.5V	Current	<500mA
Power consumption	< 0.04w	N.W.	about 400 g
Size	130*102*34 mm	Battery-operated	3×1.5v AA


### Characteristic

- Self calibration
- Works in dry or wet canals (Blood, Hypochlorite, Saline, etc)
- 0.1 mm length determination
- Colored LCD screen display with real-time graphic of root canals.
- Adjustable volume of the sound indicator

- Adjustable WL bar for different anatomies
- Steady reading around the apical foramen (apical constriction)
- Excellence in working length determination comes with 5 different setting of apical adjustment around the minor diameter
- Autoclavable file holder and contrary electrode to avoid cross infection
- Operated by 3 pieces of AA batteries to save the power

#### A. LCD Display





**B. Autoclavable accessories: file holders and contrary electrodes**

